Contour[®]TS

Jalur Ujian Glukosa Darah

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Rujukan:

- Sewell DL. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline, 3rd Edition*. Clinical and Laboratory Standards Institute. CLSI document M29-A3; ISBN 156238-567-4. March 2005.
- Burtis CA, Ashwood ER, editors. *Tietz Fundamentals of Clinical Chemistry*. 5th edition. Philadelphia, PA: WB Saunders Co; 2001;444.
- Atkin SH, et al. Fingerstick glucose determination in shock. *Annals of Internal Medicine*. 1991;114(12):1020-1024.
- McEnroe RJ, et al. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. EP7-A2, vol 25, no 27. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

 ASCENSIA Diabetes Care	 Ascensia Diabetes Care Holdings AG Peter Merian-Strasse 90 4052 Basel, Switzerland
Authorized Representative: Ethos Healthcare Sdn Bhd A-4-11 Tamarind Square Persiaran Multimedia 63000 Cyberjaya Selangor, Malaysia Tel: (60) 1 9315 9818 Email: enquiry@ethosmy.com Registration No.: IVDC1104417317	Address of Production Facility: PHC Corporation In Vitro Diagnostics Division 2131-1 Minamigata, Toon, Ehime 791-0395, Japan
www.diabetes.ascensia.com	www.patents.ascensia.com

Ascensia, logo Ascensia Diabetes Care, Contour dan logo No Coding ialah tanda dagangan dan/atau tanda dagangan berdaftar Ascensia Diabetes Care Holdings AG.

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Contour™TS

Blood Glucose Test Strips

Intended Use: CONTOUR™TS test strips are intended for self-testing by people with diabetes and by health care professionals to monitor glucose concentrations in whole blood.

Summary: The CONTOUR TS test strip is designed for use with the CONTOUR™TS blood glucose meter. The test provides a quantitative measurement of glucose in blood from 10 - 600 mg/dL (0.6 mmol/L to 33.3 mmol/L).

Storage and Handling:

- Store the strips at temperatures 9°C \rightarrow 30°C between 9°C and 30°C .
- **Store test strips in their original bottle only. Always close the lid immediately and tightly after removing a test strip.**
- **Wash and dry your hands well before handling the test strips.**

- Do not use the test strips after the  expiry date. The expiry date is printed on the bottle label and on the test strip carton.
- If the meter and/or test strips are moved from one temperature to another, allow 20 minutes for them to adjust to the new temperature before performing a blood glucose test. Your user guide will identify the appropriate operating temperature range for the CONTOUR TS meter you are using.

-  The test strips are for single use only. **Do not reuse test strips.**
- Upon opening the test strip carton ensure the strip bottle cap is securely closed. If the cap is not closed, do not use the test strips for testing. Examine the product for missing, damaged, or broken parts. Please contact your local Ascensia Diabetes Care agency for replacement parts and assistance. www.diabetes.ascensia.com

 **Test Procedure:** See your CONTOUR TS user guide and accompanying inserts for detailed testing instructions.

Test Results:

Please note the units of measure displayed by your meter. Results in mg/dL will **never** have a decimal point (e.g., $96 \frac{\text{mg}}{\text{dL}}$); results in mmol/L will **always** have a decimal point (e.g., $5.3 \frac{\text{mmol}}{\text{L}}$). If your test result is displayed in the wrong unit of measure, please contact your local Ascensia Diabetes Care agency. www.diabetes.ascensia.com

- **If your blood glucose reading is under 50 mg/dL (2.8 mmol/L) on the meter display, follow medical advice immediately.**
- **If your blood glucose reading is above 250 mg/dL (13.9 mmol/L) on the meter display, call your health care professional as soon as possible.**
- **Always consult your health care professional before changing your medication based on CONTOUR TS test results.**

Questionable or Inconsistent Results: See the CONTOUR TS user guide for problem solving. If attempts to correct a problem fail, please contact your local Ascensia Diabetes Care agency. www.diabetes.ascensia.com

Quality Control: You should perform a control test when using your meter for the first time, or when you open a new bottle or package of test strips, or if you think your meter may not be working properly, or if you have repeated unexpected blood glucose results. Only use CONTOUR™TS control solutions. These control solutions are designed specifically for use with the CONTOUR TS system. The control results should fall within the control range(s) printed on the test strip bottle label and the test strip carton. If they don't, do not use your meter for blood glucose testing until you resolve the issue.

Information For Safety

- For **[IVD]** *in vitro* diagnostic use only. External use, do not swallow.
- Potential Biohazard: Health care professionals or persons using this system on multiple patients should follow the infection control procedure approved by their facility. All products or objects which come in contact with human blood, even after cleaning, should be handled as if capable of transmitting infectious diseases. The user should follow the recommendations for the prevention of blood-borne transmissible diseases in health care settings as recommended for potentially infectious human specimens.¹
- Dispose of the test strips as medical waste or as advised by your health care professional.
- CONTOUR TS test strips contain product of animal origin BSA (Bovine Serum Albumine) as part of the non-reactive ingredients. The risk of infection from this material has been assessed to be negligible.

Chemical Composition: FAD glucose dehydrogenase (*Aspergillus sp.*, 2.0 U/test strip), 6%; potassium ferricyanide 56%; Non-reactive ingredients 38%.
Comparison Options: The CONTOUR TS system is designed for use with venous and capillary whole blood. Comparison to a laboratory method must be done simultaneously with aliquots of the same sample. Note: Glucose concentrations drop rapidly due to glycolysis (approximately 5–7% per hour).²

Limitations:

1. **Preservatives:** Blood may be collected by health care professionals into test tubes containing heparin. Do not use other anticoagulants or preservatives.
2. **Altitude:** Up to 3048 meters does not significantly affect results.
3. **Alternative Site Testing:** Not for Alternative Site Testing.
4. **Peritoneal dialysis solutions:** Icodextrin does not interfere with CONTOUR TS test strips.
5. **Xylose: Do not use during or soon after xylose absorption testing. Xylose in the blood will cause an interference.**
6. **Contraindications:** Capillary blood glucose testing may not be clinically appropriate for persons with reduced peripheral blood flow. Shock, severe hypotension, hyperosmolar hyperglycemia and severe dehydration are examples of clinical conditions that may adversely affect the measurement of glucose in peripheral blood.³
7. **Hematocrit:** CONTOUR TS test strip results are not significantly affected by hematocrit levels in the range of 0% to 70%.
8. **Interference:** The CONTOUR TS system was tested against the following potentially interfering substances occurring naturally in the blood: bilirubin, creatinine, galactose, glutathione, hemoglobin, triglycerides, and uric acid. No interfering effect was observed for these substances at the highest concentration⁴ of either the common pathological level or three times the upper reference value. **Cholesterol values of >300 mg/dL (>7.8 mmol/L) may produce results with greater than 10% bias.**
9. **Interference:** The CONTOUR TS system was tested against the following potentially interfering substances occurring from therapeutic treatments: ascorbic acid, paracetamol (acetaminophen), dopamine, sodium gentisate, ibuprofen, icodextrin, L-dopa, maltose, methyl dopa, pralidoxime iodide, sodium salicylate, tolbutamide. No interfering effect was observed for these substances at the highest concentration⁴ of either the toxic concentration or three times the maximum therapeutic concentration. **Tolazamide concentrations of >5 mg/dL (>0.16 mmol/L) may produce results with greater than 10% bias.**

The CONTOUR TS meter is not indicated for neonatal use.


References:

1. Sewell DL. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline, 3rd Edition.* Clinical and Laboratory Standards Institute. CLSI document M29-A3; ISBN 156238-567-4. March 2005.
2. Burtis CA, Ashwood ER, editors. *Tietz Fundamentals of Clinical Chemistry.* 5th edition. Philadelphia, PA: WB Saunders Co; 2001;444.
3. Atkin SH, et al. Fingerstick glucose determination in shock. *Annals of Internal Medicine.* 1991;114(12):1020-1024.
4. McEnroe RJ, et al. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition.* EP7-A2, vol 25, no 27. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

Ascensia Diabetes Care Singapore Pte Ltd
600 North Bridge Road,
#06-05 Parkview Square,
Singapore 188778
Email: contourcustomerservicesg@ascensia.com

Distributed in Singapore by:
Pharmaforte Singapore Pte Ltd
6 Tagore Drive,
#03-11, Tagore Building,
Singapore 787623
Customer service no.: (65) 6452 8488
www.diabetes.ascensia.com

Ascensia Diabetes Care Holdings AG
Peter Merian-Strasse 90
4052 Basel, Switzerland

Ascensia, the Ascensia Diabetes Care logo, Contour and the No Coding logo are trademarks and/or registered trademarks of Ascensia Diabetes Care Holdings AG.
www.patents.ascensia.com

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© 2024 Ascensia Diabetes Care Holdings AG. Hak cipta terpelihara.